PhosMap

Dongdong Zhan

2019-06-24

Introduction

PhosMap is a comprehensive R package for analyzing quantitative phosphoproteomics data, which provides mutltiple functions for users as follow: 1. clustering: principal component analysis (PCA) and t-Distributed Stochastic Neighbor Embedding (t-SNE); 2. differential expression analysis; 3. time course analysis; 4. kinase activity prediction to find activated/deactivated kinases from the kinase-substrate database 5. phosphorylation motif enrichment analysis to provide clues for finding candidate kinases that are not present in the database; 6. and data visualization.

Loading data

To test the efficacy of PhosMap and help users get started quickly, we collected a dataset including 39 phosphoproteomics and 32 proteomics raw files deposited in the ProteomeXchange Consortium (Ressa, et al., 2018). The partial key intermediate results were provided for uers to master PhosMap. We have embeded intermediate results from demo into PhosMap for help users get started.

```
# Load PhosMap
library("PhosMap")
# Load intermediate results from https://github.com/ecnuzdd/PhosMap_datasets
ftp_url <- "ftp://111.198.139.72:4000/pub/PhosMap_datasets/BRAFi.RData"
load_data <- load_data_with_ftp(ftp_url, 'Rdata')

## First loading data from FTP sever, it may take a few minutes.

## Downloading data from ftp://111.198.139.72:4000/pub/PhosMap_datasets/BRAFi.RData.

## Completing the RData load.

temp_file <- tempfile()
writeBin(load_data, temp_file)
load(temp_file)</pre>
```

Object description

- background_df A data frame for motif enrichment analysis as background
- combined df with mapped gene symbol Get a data frame mapped GI number to Gene Symbol
- data_frame_normalization_with_control_no_pair A data frame containing phosphoproteomics data normalized by proteomics data.
- foreground_df A data frame for motif enrichment analysis as foreground.
- fuzzy_input_df A data frame for time course analysis as input.
- group A factor for experiment group information.
- merge_df_with_phospho_peptides A merged phosphoproteomics data frame based on peptides files

- (unique ID).
- motif_group_m_ratio_df_mat A matrix for motif profile.
- phospho_data_filtering_STY_and_normalization A phosphoproteomics data frame after normalization and filtering.
- profiling_data_normalized A proteomics data frame after normalization and filtering.
- summary_df_of_unique_proteins_with_sites A data frame that phosphorylation sites had been mapping to protein sequence and eliminated redundancy.

Data pre-processing

An intact data pre-processing procedure of phosphoproteomics data covered three main steps: merging input files after quality control, mapping phosphorylation sites (p-sites) to the corresponding protein sequence and data normalization.

Merging input files after quality control

'Phosphoproteomics data' and 'The phosphoproteomics experimental design' are required as input. For p-sites detected by Mascot, PhosMap could provide confidence probability of p-sites extracted from Mascot xml file. For p-sites detected by other software, a two column table including their corresponding sequences and confidence probability was indispensable. Then quality control at peptide and site levels for each experiment was performed.

```
BASE_DIR <- getwd() # working directory</pre>
BASE DIR <- normalizePath(BASE DIR)
phosphorylation_exp_design_info_file_path <- normalizePath(file.path(BASE_DIR,</pre>
                                                                 'phosphorylation_exp_design_info.txt'))
phosphorylation_peptide_dir <- normalizePath(file.path(BASE_DIR, 'phosphorylation_peptide_txt'))</pre>
if(FALSE){
  # if you have xml files from mascot results, you can run the cmd to parser them to text files.
  mascot_xml_dir <- normalizePath(file.path(BASE_DIR, 'mascot_xml'))</pre>
  mascot_txt_dir <- normalizePath(file.path(BASE_DIR, 'mascot_txt'))</pre>
  extract_psites_score(phosphorylation_exp_design_info_file_path, mascot_xml_dir, mascot_txt_dir)
  # Based on above-mentioned text files from Mascot results,
  # the following cmd can generate CSV files of phosphorylation sites with confidence score.
  psites_score_dir <- normalizePath(file.path(BASE_DIR, 'psites_score_txt'))</pre>
  generate_psites_score_file(mascot_txt_dir, phosphorylation_peptide_dir, psites_score_dir)
# Merge phosphoproteomics data based on peptides files (unique ID).
# If qc = TRUE, considering confidence score of phosphorylation sites.
# A merged phosphoproteomics data frame based on peptides files (unique ID).
merge_df_with_phospho_peptides <- pre_process_filter_psites(</pre>
  phosphorylation_peptide_dir,
  psites_score_dir,
  phosphorylation_exp_design_info_file_path,
  qc = TRUE,
  min_score = 20,
  min_FDR = 0.01
)
```

Mapping phosphorylation sites (p-sites) to the corresponding protein sequence

Mapping protein gi to gene symbol and outputing expression profile matrix with gene symnol.

```
# Get a data frame mapped GI number to Gene Symbol.
# system.time({
# combinated_df_with_mapped_gene_symbol = get_combined_data_frame(
# merge_df_with_phospho_peptides, species = 'human', id_type = 'RefSeq_Protein_GI'
# )
# })
head(combined_df_with_mapped_gene_symbol)
```

```
Sequence
                                   Modification GeneSymbol Exp027012
##
## 1 AAAAAATAPPsPGPAQPGPR gi|39930517 Phospho (ST)(11) SAMD1 191101700
        AAALSLSTLAsPK gi 4758248 Phospho (ST)(11)
                                                 EFNB1 19612020
        AAALSLSTLASPK gi|4758248 Phospho (ST)(5) EFNB1 225996500

AAALSLSTLASPK gi|4758248 Phospho (ST)(7) EFNB1 25986220
## 3
## 4
## 5
        AAAtPESQEPQAK gi|13491174 Phospho (ST)(4) MARCKSL1 543292900
## 6 AADPPAENSsAPEAEQGGAE gi 34098946 Phospho (ST)(10) YBX1 385447600
  Exp027013 Exp027014 Exp027015 Exp027016 Exp027017 Exp027018
## 1 72397330 0 32192620 231611600 124170700 28650550
## 2 17807950 19357340 79808370 84246970 56323030 48652040
## 3 148918100 77889940 360054600 226245100 189362500 203163300
## 4 12633960 0
                        0
                                      0
## 5 1942228000 1524829000 1160647000 2877282000 1013506000 1597019000
## 6 624299300 579280200 364341700 649720200 415048000 890674400
  Exp027019 Exp027020 Exp027021 Exp027022 Exp027023 Exp027024
## 1 31780240 123238200 30951010 0 159176400 0
## 2 33331610
              0 26782070 59331840 81300420 11794550
## 3 159109000 112891900 117936700 130344300 136386600 40931780
    0 0 0
                                0 0
                                                       0
## 4
## 5 578654800 2225267000 1101911000 2347455000 800039800 32572150
## 6 803288700 1327916000 637003800 1049919000 917176600 96717140
##
   Exp027025 Exp027026 Exp027027 Exp027028 Exp027029 Exp027030
## 1 0 0 36968750 345538800 167434500 232864900
## 2 25698570 48053290 11260400 49733950 17026620 33251190
## 3 40848750 84373550 18245820 449307600 27093610 205768900
      0 0 0 0 0 7733698
## 4
## 5 1403377000 2010127000 1768473000 2510754000 2389508000 2335473000
## 6 408221700 880458700 551895600 1435267000 626601400 653992300
   Exp027031 Exp027032 Exp027033 Exp027034 Exp027035 Exp027036
##
          0 341993300 171803500 0 38057890 109961100
## 1
## 2 68674110 82809850 20488140 66844210 94660320 26547410
## 3 140357500 203682300 71224800 170784100 153310800 41458280
          0 5456458 0
                                     0 0
## 5 2901491000 2740149000 1173034000 315909400 341922700 1432846000
## 6 693000900 766916700 568389600 639079000 754503200 1056575000
    Exp027037 Exp027038 Exp027039 Exp027040 Exp027041 Exp027042
##
## 1 167525000 165389700 51471910 143915200 356015900 208103900
## 2 0 0 44516010 107605300 33192480 0
## 3 70873710 99639770 132810800 78469830 92153170 344426900
                       9818745 0
## 4
          0
                     0
                                                0
## 5 2104620000 2328424000 1673055000 2601683000 2413220000 1619502000
## 6 807748200 978603400 685710400 844008900 898072800 449702300
  Exp027043 Exp027044 Exp027045 Exp027046 Exp027047 Exp027048 Exp027049
                  0 54273770 63307160 83332280 162916400 120589000
## 1 101188100
      0
                  0 0 0
## 2
                                             0 22553820
                                                               0
## 3 134677700 68013420 27048770 52237850 65281570 116332000
                                                               0
```

• Constructing the data frame with unique phosphorylation site for each protein sequence.

```
# Assign psites to protein sequence.
# Unique ID: protein_gi + phosphorylation site in protein sequence.
# system.time({
# summary_df_of_unique_proteins_with_sites = get_summary_with_unique_sites(
# combinated_df_with_mapped_gene_symbol, species = 'human', fasta_type = 'refseq'
# )
# })
head(summary_df_of_unique_proteins_with_sites)
```

```
##
                 AA_in_protein AA_in_peptide
                                                  Sequence
                               s11 AAAAAATAPPsPGPAQPGPR
## gi | 39930517_s161
                     s161
## gi|4758248_s281
                       s281
                                    s5
                                             AAALsLSTLASPK
                       s283
## gi 4758248_s283
                                     s7
                                              AAALSLSTLASPK
## gi|4758248 s287
                       s287
                                   s11
                                             AAALSLSTLAsPK
                                     s3 GGsGTAGTEPSDIIIPLR
## gi 4758248_s292
                        s292
## gi 4758248_t284
                        t284
                                     t8
                                             AAALSLStLASPK
##
                        ID
                              Modification GeneSymbol Exp027012
## gi|39930517_s161 gi|39930517 Phospho (ST)(11) SAMD1 191101700
## gi|4758248_s281 gi|4758248 Phospho (ST)(5)
                                            EFNB1 225996500
## gi|4758248_s283 gi|4758248 Phospho (ST)(7)
                                            EFNB1 25986220
## gi|4758248_s287 gi|4758248 Phospho (ST)(11)
                                            EFNB1 19612020
## gi | 4758248_s292 gi | 4758248 Phospho (ST)(3)
                                            EFNB1
## gi|4758248_t284 gi|4758248 Phospho (ST)(8)
                                            EFNB1
                Exp027013 Exp027014 Exp027015 Exp027016 Exp027017
## gi|39930517_s161 72397330 0 32192620 231611600 124170700
## gi | 4758248_s281 148918100 77889940 360054600 226245100 189362500
## gi|4758248_s283 12633960 0 0
## gi|4758248_s292
                 0 0
                                     0
## gi 4758248_t284
                       0 12164570
                                        0 34897000
##
                Exp027018 Exp027019 Exp027020 Exp027021 Exp027022
## gi|39930517_s161 28650550 31780240 123238200 30951010
## gi | 4758248_s281 203163300 159109000 112891900 117936700 130344300
## gi | 4758248_s283
                  0
                           0
                                  0
                                               0
## gi|4758248 s287
                 48652040 33331610
                                        0 26782070 59331840
## gi|4758248_s292
                       0
                                0
                                        0
                                                 0
## gi 4758248_t284
                       0
                              0
                                       0
                                                 0
##
                Exp027023 Exp027024 Exp027025 Exp027026 Exp027027
## gi|39930517_s161 159176400 0 0
                                                0 36968750
## gi | 4758248_s281 136386600 40931780 40848750 84373550 18245820
## gi|4758248_s283 0
                              0
                                        0
## gi|4758248_s287 81300420 11794550 25698570 48053290 11260400
## gi | 4758248_s292
                       0
                                0
                                        0
```

```
##
              Exp027028 Exp027029 Exp027030 Exp027031 Exp027032
## gi|39930517 s161 345538800 167434500 232864900 0 341993300
## gi | 4758248 s281 449307600 27093610 205768900 140357500 203682300
## gi|4758248 s283
               0
                          0 7733698 0
                                           5456458
0
## gi|4758248_s292
                           0
                                  0
                                         0
## gi|4758248 t284 12024320
                          0
                                  0
##
              Exp027033 Exp027034 Exp027035 Exp027036 Exp027037
## gi | 4758248_s281 71224800 170784100 153310800 41458280 70873710
## gi|4758248_s283
               0 0
                               0
## gi 4758248_s287 20488140 66844210 94660320 26547410
## gi | 4758248 s292
                   0
                          0
                                  0
                                          0
## gi|4758248_t284
                   0
                          0
                                 0 6888869
              Exp027038 Exp027039 Exp027040 Exp027041 Exp027042
## gi | 39930517_s161 165389700 51471910 143915200 356015900 208103900
## gi|4758248_s281 99639770 132810800 78469830 92153170 344426900
                             0
0
## gi|4758248_s287
                  0 44516010 107605300 33192480
                  0 0 7382834 0
## gi 4758248 s292
                  0 40079320 0 35790700 21361960
## gi|4758248_t284
##
              Exp027043 Exp027044 Exp027045 Exp027046 Exp027047
## gi | 4758248_s281 134677700 68013420 27048770 52237850 65281570
                       0
                  0
## gi 4758248_s283
                                  0
                                                 0
## gi 4758248_s287
                   0
                          0
                                  0
                                          0
                                                 a
## gi 4758248_s292
                  0
                          0
                                  0
                                          0
                                                 0
## gi 4758248_t284
              8331117
                          0
                                  0
              Exp027048 Exp027049 Exp027050
## gi|39930517_s161 162916400 120589000 271443600
## gi|4758248_s281 116332000
              0
## gi|4758248_s283
                                  0
## gi | 4758248 s287 22553820
                          0
                                  0
## gi|4758248 s292
## gi 4758248_t284
```

Data normalization

PhosMap provides two kinds of normalizations.

PhosMap allowed for a total sum scaling normalization.

```
# Imputation with the next order of magnitude of the minimum except for zero.
# Filtering data only including phosphorylation site.
phospho_data_filtering_STY_and_normalization_list <- get_normalized_data_of_psites(
    summary_df_of_unique_proteins_with_sites,
    phosphorylation_exp_design_info_file_path,
    topN = NA, mod_types = c('S', 'T', 'Y')
)
phospho_data_filtering_STY <-
    phospho_data_filtering_STY_and_normalization_list$ptypes_area_df_with_id
phospho_data_filtering_STY_and_normalization_list$ptypes_fot5_df_with_id
head(phospho_data_filtering_STY_and_normalization)</pre>
```

2. If having matched proteomics data with phosphoproteomics, PhosMap allowed for normalizing phosphoproteomics data based on proteomics data.

```
# Based on phospho_data_filtering_STY_and_normalization
ID <- paste(phospho_data_filtering_STY_and_normalization$GeneSymbol,</pre>
           phospho_data_filtering_STY_and_normalization$AA_in_protein,
           sep = '_')
Value <- phospho_data_filtering_STY_and_normalization[,-seq(1,6)]</pre>
phospho_data <- data.frame(ID, Value)</pre>
phospho_data_rownames <- paste(phospho_data_filtering_STY_and_normalization$ID,</pre>
                                phospho_data_filtering_STY_and_normalization$GeneSymbol,
                                phospho_data_filtering_STY_and_normalization$AA_in_protein,
                                sep = ' ')
rownames(phospho_data) <- phospho_data_rownames</pre>
# Further normalize phosphoproteomics data based on proteomics data
# The configurations of function see help document.
data_frame_normalization_with_control_no_pair <- normalize_phos_data_to_profiling(</pre>
  phospho_data, profiling_data_normalized,
  phosphorylation_exp_design_info_file_path,
  profiling_exp_design_info_file_path,
  control_label = '0',
  pair_flag = FALSE
head(data_frame_normalization_with_control_no_pair)
```

Data analysis

PhosMap incorporated four analysis modules, including clustering and differential expression analysis, time course analysis, kinase-substrate enrichment analysis to find activated/deactivated kinases and motif enrichment analysis.

Clustering and differential expression analysis

In PhosMap, Clustering methods allowed for t-SNE and PCA


t-SNE

```
# Clustering: t-SNE or PCA
expr_data_frame <- data_frame_normalization_with_control_no_pair
# t-SNE using all experiments
visualization_with_simple_tsne(expr_data_frame, group)</pre>
```

```
## Loading required namespace: Rtsne
```

```
## Performing PCA
## Read the 39 x 39 data matrix successfully!
## OpenMP is working. 1 threads.
## Using no_dims = 2, perplexity = 10.000000, and theta = 0.500000
## Computing input similarities...
## Building tree...
## Done in 0.00 seconds (sparsity = 0.902038)!
```


```
## Learning embedding...
## Iteration 50: error is 55.494115 (50 iterations in 0.01 seconds)
## Iteration 100: error is 56.304273 (50 iterations in 0.01 seconds)
## Iteration 150: error is 60.784762 (50 iterations in 0.01 seconds)
## Iteration 200: error is 56.973069 (50 iterations in 0.00 seconds)
## Iteration 250: error is 60.232398 (50 iterations in 0.01 seconds)
## Iteration 300: error is 1.697451 (50 iterations in 0.01 seconds)
## Iteration 350: error is 1.123384 (50 iterations in 0.01 seconds)
## Iteration 400: error is 0.957553 (50 iterations in 0.00 seconds)
## Iteration 450: error is 0.807439 (50 iterations in 0.01 seconds)
## Iteration 500: error is 0.637739 (50 iterations in 0.01 seconds)
## Iteration 500: error is 0.637739 (50 iterations in 0.01 seconds)
## Iteration performed in 0.08 seconds.
```




clustering example: t-SNE

PCA

```
expr_ID <- as.vector(expr_data_frame[,1])
expr_Valule <- expr_data_frame[,-1]
expr_Valule_mean <- NULL
expr_Valule_row <- nrow(expr_Valule)
for(i in 1:expr_Valule_row){
    x <- as.vector(unlist(expr_Valule[i,]))
    x_m <- tapply(x, group, mean)
    expr_Valule_mean <- rbind(expr_Valule_mean, x_m)
}
group_levels = levels(group)
colnames(expr_Valule_mean) <- group_levels
expr_df <- data.frame(expr_ID, expr_Valule_mean)
# PCA using mean value in group for comparison with original Literature
visualization_with_simple_pca(expr_df)</pre>
```





 In PhosMap, differential expression analysis methods allowed for limma, SAM and ANOVA Data preparation (t2 VS t0)


```
# Differently expressed Proteins/Genes analysis
# t2 vs t0
expr_data_frame <- data_frame_normalization_with_control_no_pair[,1:17]
# phosphoproteomics data normalized by proteomics data
# select phosphorylation sites with greater variation
expr_data_frame_var <- apply(expr_data_frame, 1, function(x){
    var(x[-1])
})
index_of_kept <- which(expr_data_frame_var>1)
expr_data_frame <- expr_data_frame[index_of_kept,]

# group information (t0 vs t2)
deps_group_levels <- c('t0', 't2')
deps_group <- factor(as.vector(group)[1:16], levels = deps_group_levels)</pre>
```

limma

```
# (1) Limma
limma_results_df <- analysis_deps_limma(expr_data_frame, deps_group, deps_group_levels,</pre>
                                          log2_label = FALSE, adjust_method = 'none')
## Loading required namespace: limma
##
   The matrix of experiment design.
                                         groupt0 groupt2
            1
                     0
## 1
## 2
            1
                     0
## 3
            1
                     0
## 4
            1
## 5
            1
                     0
            1
                     0
## 6
## 7
            1
                     0
## 8
            1
                     0
## 9
```

```
## 10
## 11
            0
## 12
            0
                    1
            0
## 13
## 14
## 15
                    1
## 16
## attr(,"assign")
## [1] 1 1
## attr(,"contrasts")
## attr(,"contrasts")$group
## [1] "contr.treatment"
##
##
   The combination of pairwise comparison(s).
##
##
##
  The matrix of comparison statement, compare other groups with control.
                                                                                  Contrasts
##
## Levels t2-t0
##
       t0
             -1
       t2
```



Differential expression analysis: limma

SAM


ANOVA

Time course analysis

Fuzzy clustering was applied to time course analysis for discovering patterns associated with time points in PhosMap. The corresponding line chart combined with membership for each cluster was also drawn.

```
group_levels <- levels(group)
# fuzzy c-means clustering
set.seed(1000)
fuzzy_clustObj <- visualization_fuzzycluster(
    fuzzy_input_df, group, group_levels,
    k_cluster=9, iteration = 100,
    mfrow = c(3,3), min_mem = 0.1,
    plot = TRUE
)</pre>
```

```
## Loading required namespace: e1071
## Loading required namespace: ClueR
```



Time course analysis example

```
# clusters information
clusterS_info <- fuzzy_clustObj$cluster
clusterS_names <- names(clusterS_info)
clusters_df <- data.frame(clusterS_names, clusterS_info)
# write.csv(clusters_df, 'clusters_df.csv', row.names = TRUE)</pre>
```

Kinase activity prediction to find activated/deactivated kinases

In PhosMap, three kinase activity prediction methods were included: KSEA, multiple linear regression (MLR) and Mean Value.

Data preparation

```
# For early and late response
# early -> clusterS_info==1
# late -> clusterS_info==2
cluster_flag <- 'early'
cluster_symbol <- clusterS_names[clusterS_info==1]
expr_data_frame <- data_frame_normalization_with_control_no_pair</pre>
```

```
index_of_cluster <- match(cluster_symbol, expr_data_frame$ID)
cluster_df <- expr_data_frame[index_of_cluster,]</pre>
```


KSEA

```
##
## Starting KSEA
## completing: 1 / 39
## completed: 1 / 39
## completing: 2 / 39
## completed: 2 / 39
## completing: 3 / 39
## completed: 3 / 39
## completing: 4 / 39
## completed: 4 / 39
## completing: 5 / 39
## completed: 5 / 39
## completing: 6 / 39
## completed: 6 / 39
## completing: 7 / 39
## completed: 7 / 39
   completing: 8 / 39
##
## completed: 8 / 39
## completing: 9 / 39
## completed: 9 / 39
##
   completing: 10 / 39
## completed: 10 / 39
## completing: 11 / 39
## completed: 11 / 39
## completing: 12 / 39
## completed: 12 / 39
## completing: 13 / 39
## completed: 13 / 39
## completing: 14 / 39
## completed: 14 / 39
## completing: 15 / 39
## completed: 15 / 39
## completing: 16 / 39
## completed: 16 / 39
## completing: 17 / 39
## completed: 17 / 39
## completing: 18 / 39
## completed: 18 / 39
## completing: 19 / 39
## completed: 19 / 39
## completing: 20 / 39
## completed: 20 / 39
## completing: 21 / 39
## completed: 21 / 39
## completing: 22 / 39
##
   completed: 22 / 39
```

```
## completing: 23 / 39
## completed: 23 / 39
## completing: 24 / 39
## completed: 24 / 39
## completing: 25 / 39
## completed: 25 / 39
## completing: 26 / 39
## completed: 26 / 39
## completing: 27 / 39
## completed: 27 / 39
## completing: 28 / 39
## completed: 28 / 39
## completing: 29 / 39
## completed: 29 / 39
## completing: 30 / 39
## completed: 30 / 39
## completing: 31 / 39
## completed: 31 / 39
## completing: 32 / 39
## completed: 32 / 39
## completing: 33 / 39
## completed: 33 / 39
## completing: 34 / 39
## completed: 34 / 39
## completing: 35 / 39
## completed: 35 / 39
## completing: 36 / 39
## completed: 36 / 39
## completing: 37 / 39
## completed: 37 / 39
## completing: 38 / 39
## completed: 38 / 39
## completing: 39 / 39
## completed: 39 / 39
## Ending KSEA
## Extracting information data frame derived from KSEA
## ****** Regulation direction from KSEA *******
## ******* Pvalue from KSEA *******
   ****** Activity from KSEA *******
##
## ******* Kinase_site_substrate quantification matrix after KSEA ********
##
## KSEA OK! ^_^
```

```
# Activity of regulons for regulation
ksea_regulons_activity_df_cluster <- summary_df_list_from_ksea_cluster$ksea_regulons_activity_df
ksea_id_cluster <- as.vector(ksea_regulons_activity_df_cluster[,1])
ksea_value_cluster <- ksea_regulons_activity_df_cluster[,-1]
if(FALSE){
    # Pvalue of regulons for regulation
    ksea_regulons_pvalue_cluster <- summary_df_list_from_ksea_cluster$ksea_regulons_pvalue_df
    # Activity of regulons for regulation
    ksea_regulons_activity_cluster <- summary_df_list_from_ksea_cluster$ksea_regulons_activity_df
    # Expression ratio of regulons for regulation
    ksea_kinase_site_substrate_original_ratio_cluster <-
        summary_df_list_from_ksea_cluster$ksea_kinase_site_substrate_original_ratio_df
}</pre>
```

```
# plot pheatmap
if(TRUE){
  # annotation setting
  annotation_col <- data.frame(</pre>
    group = group
  )
  rownames(annotation_col) <- colnames(ksea_value_cluster)</pre>
  # breaks and colors setting
  breaks_1 <- seq(-4, -2, 0.2)
  colors_1 <- colorRampPalette(c('#11264f', '#145b7d'))(length(breaks_1)-1)</pre>
  breaks_2 <- seq(-2, -1, 0.2)
  colors_2 <- colorRampPalette(c('#145b7d', '#009ad6'))(length(breaks_2))</pre>
  breaks_3 <- seq(-1, 1, 0.2)
  colors_3 <- colorRampPalette(c('#009ad6', 'white', '#FF6600'))(length(breaks_3))</pre>
  breaks_4 <- seq(1, 2, 0.2)
  colors_4 <- colorRampPalette(c('#FF6600', 'red'))(length(breaks_4))</pre>
  breaks_5 <- seq(2, 4, 0.2)
  colors_5 <- colorRampPalette(c('red', 'firebrick'))(length(breaks_5))</pre>
  breaks <- c(breaks_1, breaks_2, breaks_3, breaks_4, breaks_5)</pre>
  breaks <- breaks[which(!duplicated(breaks))]</pre>
  color <- c(colors_1, colors_2, colors_3, colors_4, colors_5)</pre>
  color <- color[which(!duplicated(color))]</pre>
  library(pheatmap)
  ph <- pheatmap(</pre>
   ksea_value_cluster,
    scale = 'none',
    annotation_col = annotation_col,
    clustering_distance_rows = 'euclidean',
    fontsize row = 5,
    # cutree_rows = 1,
    show_rownames = TRUE,
    fontsize_col = 5,
    # cutree_cols = 1,
    cluster_cols = FALSE,
    border_color = 'black',
    cellwidth = 5, cellheight = 5,
    breaks = breaks,
    color = color,
    legend_breaks = c(-4, -2, -1, 0, 1, 2, 4),
    legend_labels = c(-4, -2, -1, 0, 1, 2, 4),
    main = paste('Kinase-substrate enrichment analysis', cluster_flag, sep = ' ')
  )
}
```



KSEA method

MLR

Mean value

Motif enrichment analysis (MEA)

PhosMap allowed for performing MEA on user defined phosphopeptides lists.

Data preparation

```
# *** foreground ***
foreground_data <- phospho_data_filtering_STY_and_normalization # pre-processed data
foreground_sequence <- as.vector(foreground_data$Sequence)
GI <- as.vector(foreground_data$GI)</pre>
```

```
Sequence <- as.vector(foreground_data$Sequence)

AA_in_protein <- as.vector(foreground_data$AA_in_protein)

# *** required parameters ***

fixed_length <- 15

species <- 'human'

motifx_pvalue <- 0.01

# get foreground data frame

# foreground_df = get_aligned_seq_for_mea(ID, Sequence, AA_in_protein, fixed_length,

# species = 'human', fasta_type = 'refseq')

# get background_df = get_global_background_df(species = 'human', fasta_type = 'refseq')
```

Motif enrichment analysis

```
# construct foreground and background
# To facilitate testing the module, select an appropriate number of items at random.
foreground <- as.vector(foreground_df$aligned_seq)
foreground <- foreground[sample(length(foreground), 1000)]
background <- as.vector(background_df$Aligned_Seq)
background <- background[sample(length(background), 10000)]

motifs_list <- mea_based_on_background(foreground, AA_in_protein, background, motifx_pvalue)</pre>
```

```
## Start executing motifx and find motif pattern.
## Foreground sequences: 1000.
## Background sequences: 10000.
## Phosphorylation: [STY] exists in foreground.
## Motifx pvalue cutoff: 0.01.
## Motifx analysis OK! ^ ^
## $S
##
                           score foreground_matches foreground_size
## 1
       ......SPS....R 37.664847
                                                  7
                                                                 820
       ....R..SPT..... 37.282931
                                                  7
## 2
                                                                 813
       ..K...SSP..... 36.564902
                                                 11
                                                                 806
## 3
## 4
       .....TSPS..... 35.837518
                                                  5
                                                                795
## 5
       ....K..SP....E. 35.615086
                                                  6
                                                                 790
     .R.....SP.....R 35.707569
                                                  5
## 6
                                                                 784
       ....P...SPT..... 35.721156
## 7
                                                  6
                                                                779
       .....SP...S.E 35.783941
## 8
                                                  6
                                                                 773
       .S.A...SP..... 34.973447
## 9
                                                  5
                                                                 767
## 10 ...E...SP.K.... 34.860592
                                                  5
                                                                 762
## 11 .....P.SP.....A 34.925900
                                                  5
                                                                 757
## 12 G..R...SP..... 34.859468
                                                  5
                                                                 752
## 13 ...LS..SP..... 34.738924
                                                  7
                                                                 747
## 14 .....P.SPT..... 34.718201
                                                  4
                                                                 740
      .....S.SP.....Q 34.837585
                                                  7
## 15
                                                                 736
## 16 A..K...SP..... 34.632052
                                                  4
                                                                 729
## 17 ..S...KSP..... 34.655472
                                                  6
                                                                 725
## 18 ....R..SP.P.... 34.560902
                                                  7
                                                                 719
## 19 ...A...SP..T... 34.378505
                                                  4
                                                                 712
      .....SP.S...R 34.418713
                                                  5
                                                                 708
## 20
## 21 P.....TSP...... 34.499510
                                                                 703
```

| ## 22RSPE 34.216657 | | | | |
|---|------------------------|-----|-----|--|
| ## 26 ESP., 34.395167 4 609 ## 25N.SP., K. 34.356699 4 4 600 ## 27 26SP., P. 34.251677 4 682 ## 27 K.ESP., P. 34.156972 6 678 ## 27 K.ESP., P. 34.165972 6 678 ## 28SP., SPT., 34.090857 3 672 ## 29SP., 16.090800 199 609 ## 31RS.S., 32.090800 199 609 ## 31RS.S., A. 48.090800 6 477 ## 32R.S., A. 48.090800 6 477 ## 32R.S., A. 48.090800 19 464 ## 31RS.S., A. 48.090800 1 464 ## 33RS.S., A. 54.66788 5 469 ## 34R.KS., 32.090800 1 464 ## 35RR.S.S., A. 64.095800 1 464 ## 35RR.S.S., A. 64.6295 4 463 ## 36L.RA.S., 36.588762 7 459 ## 37RS.S., N. 37.19222 4 452 ## 38RT.S., S., 36.291406 4 448 ## 39RS.S., N. 37.19222 4 452 ## 38RT.S., S., 36.48459 2 437 ## 40 K.R.S., S., 36.48676 3 440 ## 40 K.R.S., S., 36.48676 3 440 ## 41 VR.S.S., 36.48687 2 435 ## 42TR.S.S., 36.48687 2 435 ## 42TR.S.S., S., 36.48676 3 448 ## 43RR.S.S., 36.48687 2 435 ## 43RR.S.S., 36.48687 3 448 ## 44S.S.S., 36.48687 3 448 ## 45RS.S., 36.391406 4 488 ## 46 K.R.S.S., 38.885971 3 442 ## 47R.S.S., 38.885971 3 422 ## 48S.S.S., 38.885971 3 422 ## 49 K.R.S.S., 38.885971 3 422 ## 40 K.R.S.S., 38.885971 3 422 ## 41RR.S.S., 38.885971 3 422 ## 42RR.S.S., 38.885971 3 422 ## 44S.S.S., 38.885971 3 422 ## 45R.S.S., 38.885971 3 422 ## 46 CS.D.E. 38.885971 3 422 ## 47R.S.S., 38.885971 3 422 ## 48S.S.S., 38.885971 3 422 ## 49S.S.S., 38.885971 3 329 ## 59 R.R.S.S 39.59484 3 337 ## 51 P.S.S.S 39.59484 3 337 ## 52R.S.S 30.59484 3 329 ## 51 P.S.S.S 30.59484 3 329 ## 52R.S.S 10.639735 57 398 ## 52R.S.S 10.639735 57 398 ## 53 R.S.S.S 28.894778 3 339 ## 54M.S.S.S 28.894778 3 339 ## 55 R.S.S.S.S 28.89806 2 307 ## 55 R.S.S.S.S 28.89806 2 307 ## | ## 22RSPE 34.216657 | 4 | 700 | |
| ## 25NSPK 34.356600 4 682 ## 26ESPP 34.151677 4 682 ## 27 K.ESP 34.165872 6 678 ## 28SP 34.165872 6 678 ## 28SP 34.165872 6 669 ## 30NRS 32.080800 199 669 ## 30NRS 32.080800 2 479 ## 31RS.S.A 48080800 6 477 ## 32R.S.S.A 48080800 6 477 ## 33RSKS 45.495708 5 469 ## 33RSKS 45.495708 5 469 ## 33RSKS 45.495708 5 469 ## 37RSKS 36.168762 7 459 ## 37RS.SV. 37.12922 4 463 ## 38RT.S.S 36.291496 4 448 ## 38RS.S 36.291496 4 448 ## 39RS.S 36.59864 3 449 ## 44R.S.S 36.68874 3 449 ## 44R.S.S 36.48679 2 435 ## 44RS.S 34.48879 2 435 ## 44RS.S 34.48879 2 435 ## 44RS.S 36.52240 4 428 ## 45RS.S 38.88971 3 424 ## 47R.S.S 38.889571 3 422 ## 47R.S.S 38.62213 4 499 ## 48RS.S 49.964230 4 428 ## 48RS.S 38.622140 4 414 ## 49DS.E 33.351385 5 419 ## 47R.S.A 31.768446 5 419 ## 48RS.S 38.6090105 4 409 ## 55RS.S 30.697105 4 409 ## 57RS.S.E 33.351385 5 7 398 ## 58RS.S.E 38.1768446 5 329 ## 59RS.S.E 38.2514015 3 337 ## 57RS.S.E 38.2514015 3 337 ## 58RS.S.E 38.2514015 3 337 ## 58RS.S.E 38.2514015 3 337 ## 59RS.S.E 38.697105 4 409 ## 59RS.S.E 38.697105 4 409 ## 59RS.S.E 38.59717 3 309 ## 66RS.S.E 38.59717 3 309 ## 67RS.S.E 38.69860 2 307 ## 68 KS.S.E 38.713005 2 334 ## 57RS.E 28.192383 3 332 ## 57RS.E 28.192383 3 322 ## 68 R.R.S.D.E 28.713005 2 334 ## 57RS.S.E 28.192383 3 322 ## 68 R.R.S.D.E 28.713005 2 309 ## 68 R.R.S.S.D.E 28.713005 2 309 ## 68 R.R.S.S.D.E 28.713005 2 309 ## 78 K.R.S.S.D.E 28.713005 2 309 ## 78 R.R.S.S.D.E 28.713005 3 309 ## 78 R.R.S.S.D.E 28 | ## 23SSP.G 34.289137 | 6 | 696 | |
| ## 26ESP. P. 34.151677 | ## 24 ESPS 34.339147 | 4 | 690 | |
| ## 26ESP. P. 34.151677 | ## 25NSPK 34.356690 | 4 | 686 | |
| ## 28 | | 4 | | |
| ## 28VSPT. 34.090857 3 672 ## 29SPT. 16.000000 190 669 ## 31RS.S.A. 48.000000 2 479 ## 31RS.S.A. 48.000000 6 477 ## 31RS.S.A. 48.000000 6 477 ## 33RS.S.A. 48.000000 1 464 ## 33RS.S.A. 32.000000 1 464 ## 33RS.S.A. 32.000000 1 464 ## 33RS.S.A. 32.000000 1 464 ## 33RS.S.S. 36.446295 4 463 ## 34R.VS. 36.508762 7 459 ## 38RT.S.S.S. 36.291406 4 448 ## 37RS.S.S. 35.048764 3 440 ## 39RS.S.S. 35.068764 3 440 ## 41 VR.S.S.S. 35.084045 2 437 ## 42RS.S.S. 35.084045 2 437 ## 42RS.S.S. 35.084045 2 437 ## 43RS.S.S. 35.084045 2 437 ## 44S.S.S. 35.084045 2 437 ## 44S.S.S. 35.084045 3 440 ## 44 VRS.S.S. 35.084045 3 442 ## 45RS.S. 35.084045 3 444 ## 46 VRS.S.S. 35.084045 3 444 ## 47 VRS.S.S. 35.084045 3 444 ## 48S.S.S. 35.084045 3 444 ## 49S.S.S. 35.084045 3 444 ## 49S.S.S. 35.085052 4 442 ## 48S.S.S. 35.085052 4 442 ## 48S.S.S. 35.085052 4 442 ## 49S.S.S. 35.085052 5 433 ## 59S.S.S. 35.085052 4 444 ## 50S.S.S. 35.085052 4 444 ## 50 | | | | |
| ## 29\$P. 16.000000 190 669 ## 30 .NR.S. 32.000000 2 479 ## 32\$P. | | | | |
| ## 30 | | 3 | 672 | |
| ## 31R.S.S.A. 48.080000 6 477 ## 32R.S.S 43.080000 2 471 ## 34R.WS 32.080000 1 464 ## 34R.WS 32.080000 1 464 ## 35R.S.S 36.508762 7 459 ## 36L.RA.S 36.508762 7 459 ## 37RS.S 37.129222 4 452 ## 38RT.S 36.291496 4 448 ## 39RS.S 0, 35.958837 4 4444 ## 44 VR.S 35. 648764 3 440 ## 41 VR.S 35. 648764 3 440 ## 42TR.S 36.301496 5 2 437 ## 43RR.S 36.160552 5 433 ## 44R.S.S 37.12922 4 452 ## 44 VR. S.S 37.648679 2 435 ## 45R.S.S 37.649943 2 424 ## 47R.S.D N 32.649943 2 424 ## 47R.S N 32.649943 2 424 ## 48S.D.E 38.885971 3 422 ## 47R.S.A 31.760446 5 419 ## 49DSE.E 33.351385 5 419 ## 49DSE.E 33.351385 5 419 ## 55R.S.D 23.149045 2 341 ## 55R.S.D 23.149045 2 341 ## 55R.S.D 30.067105 4 402 ## 55R.S.D 30.067105 4 402 ## 55D.S.E 30.234177 2 339 ## 55D.S.E 30.234177 2 339 ## 57R.S.S.D 26.598274 3 327 ## 58R.S.D 32.149045 2 341 ## 59R.S.D 30.234177 2 339 ## 57R.S.D 27.229618 2 329 ## 58R.S.D 27.229618 2 329 ## 59R.S.S.D 25.635262 3 321 ## 68C.S.D 25.531374 1 322 ## 69R.S.S.D 25.538262 3 321 ## 67R.S.S.D 25.635262 3 321 ## 68C.S.D 25.784277 1 318 ## 69R.S.S.D 25.784277 1 305 ## 69R.S.S.D 25.784277 1 305 ## 69R.S.S.D 25.784279 1 305 ## 69R.S.S.D 25.784279 1 305 ## 69R.S.S.D 25.784281 2 329 ## 69R.S.S.D 25.784281 2 329 ## 69R.S.S.D 25.784281 2 329 ## 79S.D.C 25.784281 2 329 ## 79S.D.C 25.784281 2 329 ## 79S.D.C 25.784292 3 298 ## 79S.D.C 25.784292 3 298 ## 79S.S.D 25.863199 2 293 ## 75S.S.D 24.288622 3 3 | ## 29SP 16.000000 | 190 | 669 | |
| ## 32 | ## 30MRS 32.000000 | 2 | 479 | |
| ## 33 | ## 31RS.S.A 48.000000 | 6 | 477 | |
| ## 33 | ## 32RSW 32.000000 | 2 | 471 | |
| ## 34 | | 5 | | |
| ## 35RR.SF 36.446295 | | | | |
| ## 36 | | | | |
| ## 37RS.SV. 37.129222 4 452 ## 38RT.SS 36.291496 4 448 ## 39RS.SQ. 35.958837 4 444 ## 40 .K.R.S.SS. 35.958837 4 444 ## 41 V.R.S.SS. 35.648764 3 448 ## 41 V.R.S.S.S. 35.648764 3 448 ## 41 V.R.S.S.S. 35.648764 3 448 ## 42TR.S.S. 34.448679 2 435 ## 43RR.SQ. 34.160552 5 433 ## 44A.S.DE. 40.964230 4 428 ## 45RS.S 32.649943 2 424 ## 46 LSQ.E. 38.885971 3 422 ## 47R.SAA. 31.760446 5 419 ## 48GS.DE. 36.252140 4 414 ## 49DSE.E. 33.351385 5 410 ## 49DSE.E. 33.351385 5 410 ## 50 .L.R.QS 29.847178 3 405 ## 51 PSGS. 36.607165 4 402 ## 52R.S. 10.639735 57 398 ## 53 .RSQ.E. 32.149045 2 341 ## 54 .N.SQ.E. 30.234177 2 339 ## 55D.S.EE. 28.192383 3 332 ## 58 .RS.DE. 32.249045 2 341 ## 57 .RSE.E. 28.19288 3 332 ## 58 .RS.D.E. 32.259184 2 329 ## 57 .R.S.S. L. 26.590274 3 327 ## 60 K.S.S.D.S. 25.922776 2 324 ## 61 CSQ.E. 25.784727 1 318 ## 62 .R.S.S.D.E. 25.784727 1 318 ## 63 .R.S.S.E. 28.598274 3 327 ## 66 .R.S.S.D.E. 25.785651 3 312 ## 67 .H.S.S.D.E. 25.7867561 3 312 ## 67 .H.S.S.D.E. 25.786877 1 309 ## 69 T.R.S.S. 25.968986 2 307 ## 69 T.R.S.S. 25.968986 2 307 ## 67 .H.S.S.D.E. 25.786519 3 304 ## 77 .R.S.S.E. 24.863189 2 293 ## 77 .R.S.D.E. 23.179855 2 285 | | | | |
| ## 38RT.SS 36.291496 | ## 36L.RA.S 36.508762 | 7 | 459 | |
| ## 39RS.SQ. 35.958837 | ## 37RS.SV 37.129222 | 4 | 452 | |
| ## 40 | ## 38RT.SS 36.291496 | 4 | 448 | |
| ## 41 VR.S.S. 35.084045 | ## 39RS.SQ. 35.958837 | 4 | 444 | |
| ## 41 VR.S.S. 35.084045 | ## 40 .KRSS 35.648764 | 3 | 440 | |
| ## 42TR.S.S. 34.448679 2 435 ## 43RR.SQ. 34.169552 5 433 ## 44 .AS.DE. 40.964230 4 428 ## 45RS.S N 32.649943 2 424 ## 46 LSD.E. 38.885971 3 422 ## 47R.SA.A. 31.760446 5 419 ## 48GS.DE. 33.351385 5 410 ## 49DSE.E. 33.351385 5 410 ## 50 .L.R.QS. 29.847178 3 405 ## 51 PSGS. 30.607105 4 402 ## 52R.S. 10.639735 57 398 ## 53SD.E. 32.149945 2 341 ## 54SD.E. 32.149945 2 341 ## 55D.S.E. 29.154015 3 337 ## 56S.E. 29.154015 3 337 ## 57S.E. 28.192383 3 332 ## 58S.D.E. 27.229618 2 329 ## 57S.E. 28.192383 3 322 ## 60 KS.D.E. 27.229618 2 329 ## 60 KS.D.S. 25.922776 2 324 ## 61 CSD.E. 25.635262 3 321 ## 63SD.E. 25.784727 1 318 ## 64SD.E. 25.635262 3 321 ## 65S.S 26.590174 1 322 ## 67S.E. 25.635262 3 321 ## 68S.D.E. 25.786787 1 318 ## 67S.E. 25.635262 3 321 ## 68S.S.E. 25.635262 3 321 ## 69 T.R.S.S.E. 25.635262 3 321 ## 67S.E. 25.635262 3 321 ## 68S.S.E. 25.635262 3 321 ## 69 T.R.S.S.E. 25.635262 3 321 ## 67S.E.E. 24.169979 1 305 ## 67S.E.E. 25.406896 2 307 ## 67S.S.E. 25.608519 3 304 ## 77S.E.E. 24.031658 2 298 ## 77S.E.P. 23.179855 2 285 | | | | |
| ## 43RR.SQ 34.160552 | | | | |
| ## 44 | | | | |
| ## 45RS.SN 32.649943 2 424 ## 46 LSD.E 38.885971 3 422 ## 47R.SAA. 31.760446 5 419 ## 48GS.DE 36.252140 4 414 ## 49DSE.E 33.351385 5 410 ## 50R.QS 29.847178 3 405 ## 51 PSGS 30.607105 4 402 ## 52R.S 10.639735 57 398 ## 53RSD.E 32.149045 2 341 ## 54SD.E 32.149045 2 341 ## 55D.S.EE 29.154015 3 337 ## 56R.S.S.N 28.713005 2 334 ## 57R.SE.E 28.192383 3 332 ## 58SD.E 27.229618 2 329 ## 59 AGS.S 26.590274 3 327 ## 60 KS.D.S. 25.922776 2 324 ## 61 CSD.E 26.531174 1 322 ## 63SD.E 25.635262 3 321 ## 63SD.E 25.784727 1 318 ## 64SGEE 25.458843 2 317 ## 65S.S 25.639266 3 315 ## 66R.S.S.S 25.908986 2 307 ## 67S.S 25.908986 2 307 ## 67S.S 25.908986 2 307 ## 67S.S.S 25.908986 2 307 ## 67S.S.S 24.909026 3 315 ## 67S.S.S 24.909026 3 305 ## 67S.S.S 25.908986 2 307 ## 67S.S.S 24.909026 3 305 ## 67S.S.S 24.809026 3 305 ## 67S.S.S 25.908986 2 307 ## 70S.S.S 24.809026 3 305 ## 77S.S.S 24.803189 2 298 ## 77S.S.S 24.803189 2 298 ## 77S.S.S 24.803189 2 293 ## 77S.S.S 24.803189 2 293 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 46 LSD.E 38.885971 ## 47R.SAA 31.760446 ## 48GS.DE 36.252140 ## 49DSE.E 33.351385 ## 59L.R.QS 29.847178 ## 51 PSGS 30.607105 ## 52R.S 10.639735 ## 53R.SD.E 32.149045 ## 54N.SD.E 32.149045 ## 55D.S.EE 29.154015 ## 55D.S.EE 29.154015 ## 56R.S.S.N. 28.713005 2 334 ## 57R.SE.E 28.192383 ## 58S.DE 27.229618 ## 59 AGS.S.S 26.5590274 ## 60 KS.D.E 25.635262 ## 61 CSD.E 25.635262 ## 62K.LS.S 25.635262 ## 63SD.E 25.784727 ## 64SGE 25.454843 ## 65S.S.E 25.658986 ## 66R.S.S.S 26.2596896 ## 67S.S.S 24.169979 ## 77S.S.E 24.169979 ## 79 KSS.S 24.169979 ## 71 Q.S.S.S.D. 23.997708 ## 73SE.E 24.1631658 ## 73SE.E 24.4031658 ## 75SD.P.D. 23.7179855 ## 77S.SD.P.D. 23.179855 ## 77S.SD.P.D. 23.179855 ## 77S.SD.P.D. 23.179855 | ## 44AS.DE 40.964230 | 4 | 428 | |
| ## 47R.SA.A. 31.760446 5 419 ## 48GS.DE 36.252140 4 414 ## 49DSE.E 33.351385 5 410 ## 50 .L.R.QS 29.847178 3 465 ## 51 PSGS 36.607105 4 462 ## 52R.S 10.639735 57 398 ## 53 .RSD.E 32.149045 2 341 ## 54 .NSD.E 30.234177 2 339 ## 55D.S.EE 29.154015 3 337 ## 55D.S.EE 29.154015 3 337 ## 56 .R.S.S.N 28.713005 2 334 ## 57 .RSE.E 28.192383 3 332 ## 58 .RS.DE 27.229618 2 329 ## 59 AGS.S 26.590274 3 327 ## 60 KSD.E 26.553174 1 322 ## 61 CSD.E 25.635262 3 321 ## 62 .K.LS.S 25.635262 3 321 ## 64SGEE 25.784777 1 318 ## 64SGEE 25.784777 1 318 ## 65 .R.S.S 26.590266 3 315 ## 66 .R.S.S.S 25.69896 2 307 ## 69 T.R.S.S 25.668919 3 304 ## 69 T.R.S.S.S 24.169979 1 305 ## 68 .EE S.E 25.068816 2 307 ## 70 .KKSS 23.907708 3 301 ## 71 QS.S.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E 24.031658 2 298 ## 73SE.E 24.031658 2 298 ## 73SE.E 24.031658 2 298 ## 77S.S 24.128510 3 296 ## 77S.S.S 24.288622 3 285 | ## 45RS.SN 32.649943 | 2 | 424 | |
| ## 48GS.DE 36.252140 | ## 46 LSD.E 38.885971 | 3 | 422 | |
| ## 49DSE.E. 33.351385 | ## 47RSAA 31.760446 | 5 | 419 | |
| ## 49DSE.E. 33.351385 | ## 48GS.DE 36.252140 | 4 | 414 | |
| ## 50 | | 5 | | |
| ## 51 PSGS | | | | |
| ## 52 | | | | |
| ## 53 | | | | |
| ## 54 .NSD.E 30.234177 | | 57 | 398 | |
| ## 55D.S.EE 29.154015 3 337 ## 56R.S.S.N 28.713005 2 334 ## 57RSE.E 28.192383 3 332 ## 58S.DE 27.229618 2 329 ## 60 KS.D 26.590274 3 327 ## 61 CSD.E 26.531174 1 322 ## 62K.LS.S 25.635262 3 321 ## 63SD.E 25.784727 1 318 ## 64SGEE 25.454843 2 317 ## 65SS.S G. 24.900026 3 315 ## 66R.S.SP 25.576561 3 312 ## 67H.S.D.N. 24.291578 2 309 ## 68EES.E 25.068986 2 307 ## 69 T.R.S.S 24.169979 1 305 ## 70 KKS.S 25.068519 3 304 ## 71 QS.S.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.S.S.S 24.1863189 2 293 ## 75SDE.P. 23.179855 2 285 | ## 53RSD.E 32.149045 | 2 | 341 | |
| ## 56 | ## 54NSD.E 30.234177 | 2 | 339 | |
| ## 57 | ## 55D.S.EE 29.154015 | 3 | 337 | |
| ## 58 | ## 56RS.SN 28.713005 | 2 | 334 | |
| ## 58 | ## 57RSE.E 28.192383 | 3 | | |
| ## 59 AGS.S 26.590274 ## 60 KS.DS. 25.922776 ## 61 CSD.E 26.531174 ## 62K.LS.S 25.635262 ## 63 .RSD.E 25.784727 ## 64SGEE 25.454843 ## 65 .RS.S G. 24.900026 ## 66R.S.S.P 25.576561 ## 67HS.D.N. 24.291578 ## 68EES.E 25.068986 ## 69 T.R.S.S 24.169979 ## 69 T.R.S.S 24.169979 ## 70 .KKS.S 25.068519 ## 71 QS.S.D 23.907708 ## 72D.ES.E 24.031658 ## 73SE.E.Q. 24.128510 ## 74 YR.SS 24.863189 ## 75SDS.Q. 23.671932 ## 76ESSA 24.288622 ## 77S.DE.P. 23.179855 ## 282 328 ## 77S.DE.P. 23.179855 ## 283 288 ## 77S.DE.P. 23.179855 | | | | |
| ## 60 KS.DS. 25.922776 | | | | |
| ## 61 CSD.E 26.531174 | | | | |
| ## 62K.LS.S 25.635262 3 321 ## 63 .RSD.E 25.784727 1 318 ## 64SGEE 25.454843 2 317 ## 65 .RS.SG. 24.900026 3 315 ## 66R.S.SP 25.576561 3 312 ## 67S.D.N. 24.291578 2 309 ## 68EES.E 25.068986 2 307 ## 69 T.R.S.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QS.S.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.S 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 63 | | 1 | 322 | |
| ## 64SGEE 25.454843 2 317 ## 65 .RS.SG. 24.900026 3 315 ## 66R.S.SP 25.576561 3 312 ## 67S.D.N. 24.291578 2 309 ## 68EESE 25.068986 2 307 ## 69 T.RS.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QS.S.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | ## 62K.LS.S 25.635262 | 3 | 321 | |
| ## 65 | ## 63 .RSD.E 25.784727 | 1 | 318 | |
| ## 65 | ## 64SGEE 25.454843 | 2 | 317 | |
| ## 66R.S.SP 25.576561 3 312 ## 67HS.D.N. 24.291578 2 309 ## 68EES.E 25.068986 2 307 ## 69 T.R.S.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QS.S.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | 3 | 315 | |
| ## 67HS.DN. 24.291578 2 309 ## 68EES.E 25.068986 2 307 ## 69 T.RS.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QSS.D 23.907708 3 301 ## 72DES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 68EESE 25.068986 2 307 ## 69 T.RS.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QSS.D 23.907708 3 301 ## 72DES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 69 T.RS.S 24.169979 1 305 ## 70 .KKS.S 25.068519 3 304 ## 71 QSS.D 23.907708 3 301 ## 72DES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDSQ. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 70 .KKS.S 25.068519 3 304 ## 71 QSS.D 23.907708 3 301 ## 72D.ES.E 24.031658 2 298 ## 73SE.E.Q 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDSQ 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P 23.179855 2 285 | | | | |
| ## 71 QSS.D 23.907708 3 301 ## 72DES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | | | |
| ## 72DES.E 24.031658 2 298 ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDS.Q. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | 3 | 304 | |
| ## 73SE.E.Q. 24.128510 3 296 ## 74 YR.SS 24.863189 2 293 ## 75SDSQ. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | ## 71 QSS.D 23.907708 | 3 | 301 | |
| ## 74 YR.SS 24.863189 2 293 ## 75SDSQ. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | ## 72DES.E 24.031658 | 2 | 298 | |
| ## 74 YR.SS 24.863189 2 293 ## 75SDSQ. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | ## 73SE.E.Q 24.128510 | 3 | 296 | |
| ## 75SDSQ. 23.671932 3 291 ## 76ESSA 24.288622 3 288 ## 77S.DE.P. 23.179855 2 285 | | 2 | | |
| ## 76ESSA 24.288622 3 288 ## 77S.DE.P 23.179855 2 285 | | | | |
| ## 77S.DE.P. 23.179855 2 285 | | | | |
| | | | | |
| ## 78S.DSE 24.055972 3 283 | | | | |
| | ## 78S.DSE 24.055972 | 3 | 283 | |
| | | | | |

```
##
  79
       ..V....S.S..D.. 23.041296
                                                   2
                                                                  280
##
  80
       ......S..G.... 4.616931
                                                  36
                                                                  278
       .S.....S.ED.... 23.555099
                                                   3
##
  81
                                                                  242
       R....S.S.....I. 23.592586
##
  82
                                                   1
                                                                  239
       K....S.S.....P 23.162563
                                                   3
##
  83
                                                                  238
                                                   3
       ....K..S.E...E. 22.457194
##
  84
                                                                  235
##
  85
       ...AS..S.V..... 23.264145
                                                   2
                                                                  232
       .....S.EE..A. 22.542927
                                                   2
##
  86
                                                                  230
##
  87
       ....Q.SS..S.... 23.200258
                                                   3
                                                                  228
       ..R.A..S.S..... 22.603529
                                                   3
##
  88
                                                                  225
                                                   3
##
       ...RS..SD..... 22.515248
  89
                                                                  222
       ...AS..S..E.... 22.098104
                                                   2
##
  90
                                                                  219
##
  91
       H....S.S.....I. 21.930574
                                                   1
                                                                  217
##
  92
       .....EDS.D..... 21.744894
                                                   3
                                                                  216
##
  93
       .N...S.S..F.... 21.582617
                                                   2
                                                                  213
       ..M....SD.E.... 21.782518
                                                   1
##
  94
                                                                  211
       ...H...SL.S.... 21.401080
                                                   2
##
  95
                                                                  210
## 96
      ......S..E.... 3.277582
                                                  32
                                                                  208
## 97
       .....S..S.... 3.546496
                                                  40
                                                                  176
  98
      ...K.W.S..... 19.115805
                                                   1
                                                                  136
       ....K..S.....W 18.844808
## 99
                                                   1
                                                                  135
## 100 .R.K...S.....P 22.197498
                                                   2
                                                                  134
## 101 QS....DS...... 20.950690
                                                   2
                                                                  132
## 102 .KS....S.S..... 21.369730
                                                   3
                                                                  130
## 103 ....K..SLT..... 21.335552
                                                   2
                                                                  127
## 104 QR.S...S..... 21.099034
                                                   1
                                                                  125
## 105 K.HK...S..... 20.811626
                                                   1
                                                                  124
## 106 K....S.S....I. 20.843039
                                                   1
                                                                  123
## 107 ....KR.S....S. 21.740975
                                                   2
                                                                  122
## 108 .V....GS..D.... 20.948026
                                                   2
                                                                  120
                                                   2
## 109 ....Q..SSP..... 21.937718
                                                                  118
                                                   2
## 110 R...S..S.P..... 20.896746
                                                                  116
## 111 ...S...S.A..S.. 20.640245
                                                   2
                                                                  114
## 112 ...... 2.345906
                                                  21
                                                                  112
## 113 ......K. 2.271358
                                                  12
                                                                   91
## 114 ......S.DI...I 21.039996
                                                   2
                                                                   79
## 115 E.....SSP..... 21.349231
                                                   1
                                                                   77
## 116 ....K..S..M.... 18.387706
                                                   1
                                                                   76
## 117 ...S...AG. 20.583493
                                                   2
                                                                   75
## 118 ...R...SS.....S 21.074320
                                                   2
                                                                   73
## 119 .T....ES.....S. 20.433519
                                                   2
                                                                   71
## 120 .....S...S... 2.149138
                                                  14
                                                                   69
## 121 ...... 2.051214
                                                   9
                                                                   55
## 122 .....S.S...... 2.566943
                                                  11
                                                                   46
  123 ......SDN..... 18.155014
                                                                   35
       background_matches background_size fold_increase
##
## 1
                        0
                                      5092
                                                     Inf
## 2
                        0
                                      5092
                                                     Inf
## 3
                        0
                                      5092
                                                     Inf
## 4
                        0
                                      5092
                                                     Inf
## 5
                        0
                                      5092
                                                     Inf
## 6
                        0
                                      5092
                                                     Inf
## 7
                        0
                                      5092
                                                     Inf
## 8
                        0
                                                     Inf
                                      5092
## 9
                        0
                                      5092
                                                     Inf
## 10
                        0
                                      5092
                                                     Inf
## 11
                        0
                                      5092
                                                     Inf
```

| | 12 | | Inf |
|----|----|----------|-----------------|
| | 13 | | Inf |
| | 14 | | Inf |
| | 15 | | Inf |
| | 16 | | Inf |
| | 17 | | Inf |
| | 18 | | Inf |
| | 19 | | Inf |
| | 20 | | Inf |
| | 21 | | Inf |
| | 22 | | Inf |
| | 23 | | Inf |
| | 24 | | Inf |
| | 25 | | Inf |
| | 26 | | Inf |
| | 27 | | Inf |
| | 28 | | Inf |
| ## | 29 | 336 5092 | 4.304043 |
| ## | 30 | 0 4756 | Inf |
| | 31 | | Inf |
| ## | 32 | 0 4756 | Inf |
| ## | 33 | 8 0 4756 | Inf |
| ## | 34 | 0 4756 | Inf |
| ## | 35 | 0 4756 | Inf |
| ## | 36 | | Inf |
| | 37 | | Inf |
| | 38 | | Inf |
| | 39 | | Inf |
| | 40 | | Inf |
| | 41 | | Inf |
| | 42 | | Inf |
| | 43 | | Inf |
| | 44 | | Inf |
| | 45 | | Inf |
| | 46 | | Inf |
| | 47 | | Inf |
| | 48 | | Inf |
| | 49 | | Inf |
| | 50 | | Inf |
| | 51 | | Inf |
| | 52 | | 2.681637 |
| | 53 | | 2.001037 Inf |
| | 54 | | Inf |
| | 55 | | Inf |
| | 56 | | Inf |
| | 57 | | Inf |
| | | | |
| | 58 | | Inf |
| | 59 | | Inf |
| | 60 | | Inf |
| | 61 | | Inf |
| | 62 | | Inf |
| | 63 | | Inf |
| | 64 | | Inf |
| | 65 | | Inf |
| | 66 | | Inf |
| | 67 | | Inf |
| ## | 68 | 8 0 4502 | Inf |
| | | | |
| | | | |

| 0 0 0 0 0 0 0 0 0 276 0 0 0 | 4502 4502 4502 4502 4502 4502 4502 4502 | Inf |
|---|--|--|
| 0 0 0 0 0 0 0 0 276 0 0 0 0 | 4502 4502 4502 4502 4502 4502 4502 4502 | Inf |
| 0 0 0 0 0 0 0 276 0 0 0 | 4502 4502 4502 4502 4502 4502 4502 4502 | Inf |
| 0 0 0 0 0 0 276 0 0 0 0 | 4502 4502 4502 4502 4502 4502 4502 4502 | Inf |
| 0 0 0 0 0 276 0 0 0 0 | 4502 4502 4502 4502 4502 4502 4502 4226 4226 4226 4226 4226 | Inf |
| 0 0 0 0 276 0 0 0 0 | 4502 4502 4502 4502 4502 4502 4226 4226 4226 4226 4226 | Inf |
| 0 0 0 276 0 0 0 0 0 | 4502 4502 4502 4502 4502 4226 4226 4226 4226 4226 | Inf Inf Inf Inf 2.112293 Inf Inf Inf |
| 0 0 276 0 0 0 0 0 | 4502 4502 4502 4502 4226 4226 4226 4226 4226 | Inf Inf Inf 2.112293 Inf Inf Inf |
| 0 0 276 0 0 0 0 0 | 4502 4502 4502 4226 4226 4226 4226 4226 | Inf Inf 2.112293 Inf Inf Inf Inf |
| 0 276 0 0 0 0 0 0 | 4502 4502 4226 4226 4226 4226 4226 | Inf 2.112293 Inf Inf Inf Inf |
| 276 0 0 0 0 0 0 0 | 4502 4226 4226 4226 4226 4226 | 2.112293 Inf Inf Inf Inf |
| 0 0 0 0 0 0 | 4226 4226 4226 4226 4226 | Inf Inf Inf Inf |
| 0 0 0 0 0 | 4226 4226 4226 4226 | Inf Inf Inf |
| 0 0 0 0 0 | 4226 4226 4226 | Inf Inf |
| 0 0 0 0 | 4226 4226 | Inf |
| 0 0 0 | 4226 4226 | Inf |
| 0 0 0 | 4226 | |
| 0 0 0 | | T111 |
| 0 0 | | Inf |
| 0 | 4226 | Inf |
| | 4226 | Inf |
| A | | Inf |
| | | Int Inf |
| | | |
| | | Inf |
| | | Inf |
| | | 1.868258 |
| | | 1.759209 |
| | | Inf |
| 0 | | Inf |
| 0 | 3377 | Inf |
| | | 1.846028 |
| | | 2.367384 |
| | | Inf |
| | | |
| | | Inf |
| | | Inf |
| | | 2.076087 |
| | | 2.563636 |
| | | 2.583463 |
| 0 | 2196 | Inf |
| | | |
| | | |
| | 0 0 0 0 0 0 0 0 343 169 0 0 0 0 | 0 4226 0 4226 0 4226 0 4226 0 4226 348 4226 501 3878 0 3377 0 365 0 2865 0 2865 0 2865 0 2865 |

```
motif
                          score foreground_matches foreground_size
##
     D....P.TP..... 36.676650
                                                  3
  1
                                                                 178
     ..A.G..TP..... 36.251038
                                                  4
                                                                 175
      ..RA...TP..... 35.905342
##
                                                  3
                                                                 171
     E....A.TP..... 35.896627
                                                  3
                                                                 168
     .....TPP..R.. 36.099525
                                                  4
##
  5
                                                                 165
##
  6
     ..V....TP.K.... 36.116506
                                                  3
                                                                 161
     ....KP.TP..... 34.921233
                                                  3
                                                                 158
  8 A.....TP...N.. 34.437500
                                                  2
                                                                 155
  9 S.....TP....R. 34.416099
                                                  4
                                                                 153
                                                  5
  10 .....EKTP...... 34.615805
                                                                 149
  11 E.....TP.K.... 31.621434
                                                  2
                                                                 144
## 12 .....TP..SR.. 30.375660
                                                  3
                                                                 142
## 13 ..R..S.T.D.... 32.075889
                                                  3
                                                                 139
## 14 ....G..TPP..... 28.623742
                                                  3
                                                                 136
## 15 A....SDT..... 29.743913
                                                  3
                                                                 133
## 16 ...... 9.258404
                                                 33
                                                                 130
## 17 ..R..SAT..... 27.376021
                                                  2
                                                                  97
## 18 ...R.S.T.....E. 25.158896
                                                  2
                                                                  95
## 19 IS...S.T..... 23.941215
                                                  2
                                                                  93
## 20 ...N.S.T.P..... 23.075485
                                                  2
                                                                  91
## 21 ...R.S.T.S..... 23.052091
                                                  2
                                                                  89
## 22 .R.SS..T..... 23.955142
                                                  2
                                                                  87
## 23 G..Q...T..S.... 23.835249
                                                  2
                                                                  85
## 24 .....ST..L...D 23.364351
                                                  3
                                                                  83
## 25 R...S..T.S..... 22.590859
                                                  2
                                                                  80
## 26 .D....T.P....R 22.392211
                                                  2
                                                                  78
## 27 .....SDT....G.. 21.769318
                                                  2
                                                                  76
## 28 .....ST..SE... 22.902286
                                                  3
                                                                  74
## 29 .....TE.E.D.. 22.121802
                                                  2
                                                                  71
                                                  3
## 30 .....T.SP...K 21.798585
                                                                  69
## 31 ..SI...T...S... 22.174881
                                                  2
                                                                  66
## 32 ..EK...T.P..... 21.653617
                                                  2
                                                                  64
## 33 ....SP.T.W.... 22.325458
                                                  2
                                                                  62
## 34 .....T.S.SD.. 20.945402
                                                  2
                                                                  60
## 35 ...H..GT..... 18.484326
                                                  1
                                                                  58
## 36 ......ET.....DG 20.833959
                                                  2
                                                                  57
## 37 ..I...DT..E.... 21.627701
                                                  2
                                                                  55
## 38 A.....GTA...... 20.710249
                                                  1
                                                                  53
## 39 .....S.T..... 2.176177
                                                 11
                                                                  52
## 40 ....S.ST....M. 20.995448
                                                  2
                                                                  41
## 41 A.....ST....D.. 20.315962
                                                  1
                                                                  39
## 42 .....T.M..P.. 18.050606
                                                  1
                                                                  38
## 43 .N.F.K.T..... 20.790071
                                                  2
                                                                  37
  44 ...SS.GT..... 21.044632
                                                                  35
      background_matches background_size fold_increase
##
## 1
                       0
                                     3280
                                                    Inf
## 2
                       0
                                                    Inf
                                     3280
## 3
                       0
                                                    Tnf
                                     3280
## 4
                       0
                                     3280
                                                    Inf
## 5
                       0
                                     3280
                                                    Inf
## 6
                       0
                                     3280
                                                    Inf
## 7
                       0
                                     3280
                                                    Inf
                       0
                                                    Inf
## 8
                                     3280
## 9
                       0
                                     3280
                                                    Inf
                       0
## 10
                                     3280
                                                    Inf
## 11
                       0
                                     3280
                                                    Inf
```

```
## 12
                                                    Inf
                       0
                                     3280
## 13
                       0
                                     3280
                                                    Inf
## 14
                                     3280
                                                    Inf
## 15
                       0
                                     3280
                                                    Inf
                     249
                                              3.343837
## 16
                                     3280
                       0
                                     3031
                                                    Inf
## 17
## 18
                       0
                                     3031
                                                    Inf
## 19
                                     3031
                                                    Inf
## 20
                       0
                                     3031
                                                    Inf
                       0
## 21
                                     3031
                                                    Inf
                       0
## 22
                                     3031
                                                    Inf
## 23
                       0
                                     3031
                                                    Inf
## 24
                       0
                                     3031
                                                    Inf
## 25
                       0
                                     3031
                                                    Inf
## 26
                       0
                                     3031
                                                   Inf
                                     3031
## 27
                       0
                                                   Inf
## 28
                       0
                                     3031
                                                    Inf
## 29
                       0
                                     3031
                                                   Inf
## 30
                       0
                                     3031
                                                    Inf
## 31
                       0
                                     3031
                                                   Inf
                       0
## 32
                                     3031
                                                   Inf
## 33
                       0
                                     3031
                                                    Tnf
## 34
                       0
                                     3031
                                                   Tnf
                       0
## 35
                                     3031
                                                    Inf
## 36
                       0
                                     3031
                                                   Inf
                       0
## 37
                                     3031
                                                    Inf
## 38
                       0
                                     3031
                                                    Tnf
## 39
                     277
                                     3031
                                             2.314704
## 40
                       0
                                     2754
## 41
                                     2754
                                                   Inf
## 42
                       0
                                     2754
                                                    Inf
                       0
## 43
                                     2754
                                                    Inf
## 44
                       0
                                     2754
                                                    Tnf
# Find sequences in foreground that are mapped to specific motif
foreground_sequences_mapped_to_motifs <- get_foreground_seq_to_motifs(motifs_list, foreground)</pre>
##
```


Plot motif logo

```
# The data can be used for ploting logo of sepcific motif: foreground_sequences_mapped_to_motifs
# ploting logo: Q.....S......
require(ggseqlogo)
```

```
## Loading required package: ggseqlogo

## Registered S3 methods overwritten by 'ggplot2':
## method from
## [.quosures rlang
## c.quosures rlang
## print.quosures rlang
```

```
ggseqlogo(foreground_sequences_mapped_to_motifs[[1]])
```



Plot motif logo


```
if(TRUE){
    # batch plot and count peptides for each motif
    foreground_sequences_mapped_to_motifs_count <- length(foreground_sequences_mapped_to_motifs)
    motifs <- names(foreground_sequences_mapped_to_motifs)
    peptides_count <- NULL
    for(i in seq_len(foreground_sequences_mapped_to_motifs_count)){
        l_i <- foreground_sequences_mapped_to_motifs[[i]]
        peptides_count <- c(peptides_count, length(l_i))
}</pre>
```

```
motifs_peptides_count_df <- data.frame(motifs, peptides_count)
# quantile(peptides_count, seq(0,1,0.05))
if(FALSE){
   plot_seqlogo(BASE_DIR, foreground_sequences_mapped_to_motifs, plot_min_seqs = 25)
}</pre>
```

Assign quantitative values of peptides to their motif

```
# Select motifs at least having 50 peptides
# Assign quantitative values of peptides to their motif
foreground_value <- foreground_data[,-c(seq(1,6))]</pre>
min_seqs <- 50
index_of_motifs <- which(peptides_count>=min_seqs)
motif_group_m_ratio_df <- NULL</pre>
for(i in index_of_motifs){
  motif <- motifs[i]</pre>
  aligned_peptides <- foreground_sequences_mapped_to_motifs[[i]]</pre>
  index_of_match <- match(aligned_peptides, foreground_df$aligned_seq)</pre>
  motif_value <- foreground_value[index_of_match,]</pre>
  motif_value_colsum <- colSums(motif_value)</pre>
  motif_group_m <- tapply(motif_value_colsum, group, mean)</pre>
  motif_group_m_ratio <- motif_group_m/motif_group_m[1]</pre>
  motif_group_m_ratio_df <- rbind(motif_group_m_ratio_df, motif_group_m_ratio)</pre>
}
motifs_subset <- motifs[index_of_motifs]</pre>
peptides_count_subset <- peptides_count[index_of_motifs]</pre>
rownames(motif_group_m_ratio_df) <- paste(motifs_subset, peptides_count_subset)</pre>
# The matrix is import inot pheatmap
motif_group_m_ratio_df_mat <- as.matrix(motif_group_m_ratio_df)</pre>
# plot pheatmap
if(TRUE){
  library(pheatmap)
  # breaks and colors setting
  breaks_1 <- seq(0, 0.5, 0.1)
  colors_1 <- colorRampPalette(c('green', 'blue'))(length(breaks_1)-1)</pre>
  breaks_3 <- seq(0.5, 1.5, 0.1)
  colors_3 <- colorRampPalette(c('blue', 'white', '#FFBFBF'))(length(breaks_3))</pre>
  breaks_4 <- seq(1.5, 2, 0.1)
  colors_4 <- colorRampPalette(c('#FFBFBF', 'red'))(length(breaks_4))</pre>
  breaks_5 <- seq(2, 4, 0.1)
  colors_5 <- colorRampPalette(c('red', 'firebrick'))(length(breaks_5))</pre>
  breaks <- c(breaks_1, breaks_3, breaks_4, breaks_5)</pre>
  breaks <- breaks[which(!duplicated(breaks))]</pre>
  colors <- c(colors_1, colors_3, colors_4, colors_5)</pre>
  colors <- colors[which(!duplicated(colors))]</pre>
  length(breaks)
```

```
length(which(!duplicated(colors)))
  ph <- pheatmap(</pre>
    motif_group_m_ratio_df_mat,
    scale = 'none',
    # annotation_col = annotation_col,
    clustering_distance_cols = 'euclidean',
    fontsize_row = 6, cutree_rows = 1, show_rownames = TRUE, cluster_rows = TRUE,
    fontsize_col = 6, cutree_cols = 1, show_colnames = TRUE, cluster_cols = FALSE,
    border_color = 'black',
    # color = colors,
    cellwidth = 12, cellheight = 12,
    breaks = breaks,
    color = colors,
    legend_breaks = c(0, 0.5, 1, 1.5, 2, 4),
    legend_labels = c(0, 0.5, 1, 1.5, 2, 4),
    main = 'Motif enrichment analysis'
  )
}
```



KSEA method

Formatting output

```
# formatting output
formatted_output_df <- formatted_output_mef_results(foreground_sequences_mapped_to_motifs)
## Output formatted sequences in foreground that are mapped to specific motif.</pre>
```

```
# write file
# write.table(formatted_output_df, 'formatted_output_df.txt', row.names = FALSE,

# col.names = FALSE, sep = '\t')
```

Session Info

sessionInfo()

```
## R version 3.6.0 (2019-04-26)
## Platform: x86_64-w64-mingw32/x64 (64-bit)
## Running under: Windows 10 x64 (build 17134)
## Matrix products: default
##
## locale:
## [1] LC_COLLATE=Chinese (Simplified)_China.936
## [2] LC_CTYPE=Chinese (Simplified)_China.936
## [3] LC_MONETARY=Chinese (Simplified)_China.936
## [4] LC_NUMERIC=C
## [5] LC_TIME=Chinese (Simplified)_China.936
##
## attached base packages:
## [1] stats
              graphics grDevices utils
                                          datasets methods
                                                               base
##
## other attached packages:
## [1] ggseqlogo_0.1
                     pheatmap_1.0.12 PhosMap_0.99.33
## loaded via a namespace (and not attached):
## [1] Rcpp_1.0.1 compiler_3.6.0 pillar_1.4.1
## [4] RColorBrewer_1.1-2 highr_0.8
                                           plyr_1.8.4
## [7] bitops_1.0-6 class_7.3-15
                                          tools_3.6.0
## [10] digest_0.6.19
                       evaluate_0.14
                                           Rtsne_0.15
## [13] tibble_2.1.3
                         gtable_0.3.0
                                           pkgconfig_2.0.2
## [16] rlang_0.3.4
                       yaml_2.2.0
                                           parallel_3.6.0
                         e1071_1.7-2
## [19] xfun_0.7
                                           dplyr_0.8.1
                       knitr_1.23
## [22] stringr_1.4.0
                                          tidyselect_0.2.5
                                           impute_1.58.0
## [25] grid_3.6.0
                       glue_1.3.1
## [28] R6_2.4.0
                         rmarkdown_1.13
                                           limma_3.40.2
## [31] purrr_0.3.2
                         ggplot2_3.1.1
                                           ClueR 1.4
## [34] magrittr_1.5 scales_1.0.0
                                           htmltools_0.3.6
## [37] assertthat_0.2.1 colorspace_1.4-1
                                           labeling_0.3
## [40] stringi_1.4.3
                       RCurl_1.95-4.12
                                           lazyeval_0.2.2
## [43] munsell_0.5.0
                         crayon_1.3.4
```